

Standard Specifications

OB

The OB Series' molded headrests are ergonomically formed and are valued for their robust performance at an economical price.

The OB Headrest is constructed with a metal shell layered with closed cell foam and coated with a durable rubber casing. Perfect as a transport headrest, a temporary solution, or individuals who rarely rely on a headrest.



STL-OTLG150



STL-OTSM250



QUICK REFERENCE • HARDWARE

Headrest mounts to hardware with a 1" Ball

Combination part numbers for the OB with Hardware



	SIZE	PAD ONLY	HMO400	HMO475	TWB480	TWB485
Headrest	LARGE 10.5" Long x 6" High	<u>OTLG150</u>	<u>STL-OTLG100</u>	<u>STL-OTLG175</u>	<u>STL-OTLG180</u>	<u>STL-OTLG185</u>
	SMALL 9.25" Long x 3.5" High	<u>OTSM250</u>	<u>STL-OTSM200</u>	<u>STL-OTSM275</u>	<u>STL-OTLG280</u>	<u>STL-OTSM285</u>

OTTOS9 Swing-Away



As Shown:
Swing-Away ([OTTOS9-B] Right & Left) with a Long Lateral Pad (930) on the right and a Green Egg Switch (ESGRN) on the left.

Swing-Away arms offer a sleek, unobtrusive design while featuring an easy, one-hand release mechanism for swinging pads and switches back a full 100 degrees. This ensures convenience and is important for transfers, therapy and fatigue issues. Swing-Away pads or switches can be quickly and easily locked back into their original setting.

F8675
(Figure 8)



Includes a standard 9" rod and egg switch mount.

Figure 8 mounts are uniquely designed to provide 240 degree rotation capability for a second pad, switch or embedded speaker on a traditional Swing-Away.

OTTOS9-B

Includes Standard 9" rods and choice of pads or Egg Switch Mounts



Facial Lateral Rod	
SU Rods	
Model	Size
SU260	Standard 9"
SU262	Short 6"
SU264	13"
SU265	18"

To Order Call:
1-800-965-9229
and for more info log on to:
www.stealthproducts.com

Facial Lateral Options

910 Pad	900 Pad	920 Pad (Bendable)	930 Pad (Bendable)	ESM Egg Switch Mount
3" L x 2" H	4.5" L x 2" H	6.375" L x 2.25" H	6.375" L x 1.5" H	1.75" L x 1.25" H


Swing-Away Options

<p>SM600 Mini proportional joystick mount</p>	<p>SM630 ASL MEC joystick mount</p>	<p>FO660 20" Flex fiber optic mount - Cable Included -</p>	<p>MC210 Tash min-cup switch</p>	<p>ML210 Micro Lite switch</p>
<p>FO650 Lateral rod fiber optic mount - Cable Included -</p>	<p>CN600 12" Flex micro-lite switch mount</p>	<p>SP625 Sip & Puff mount</p>	<p>EGG SWITCHES</p> <p>ESYEL ESRED ESBLK</p> <p>ESGRN ESBLU</p> <p>Also available: ES-EXT (Egg switch extension cable)</p>	

Unique Configurations

For pricing, refer to underlined reference numbers in the COMBO Section of the Pricelist.

Creative Configurations For Simple Solutions



OTTOS9-B Both sides
OTTOS9-R Right side only
OTTOS9-L Left side only
(See Swing-away page for additional ideas)

Swing-Aways

Facial lateral Swing-Aways for pads, switches, speaker sites and chin control devices can be added to the Combo.


Pad Sizes: • 900 • 910 • 920 • 930	Type of Foam/Gel: • Open Cell Foam • Closed Cell Foam • Action Gel	Type of Cover: • Dartex • Reverse Dartex • Vinyl • Sheep Skin • Neoprene
--	---	---

(Zippers available upon request)

CALL AND PLACE YOUR ORDER TODAY!
1.800.965.9229

OB01

Mini Proportional Joystic Mount



Includes: Small OB Headrest (**STL-OTS250**) with Multi Axis Hardware (**HMO475**), Swing-Away ([**OTTOS9-L**] Left) and Mini Proportional Joystick Mount (**SM600**).

(Joystick not included)

OB02

Figure 8 Bracket



Includes: Small OB Headrest (**STL-OTS250**) with Link Hardware (**TWB480-1"**) and Flip Down Mount (**FDM380**), Swing-Away ([**OTTOS9-L**] Left) with Figure 8 Bracket (**F8675**), Long Facial Lateral Pad (**920**) and Yellow Egg Switch (**ESYEL**).

OB03

Fiber Optic Mount



Includes: Small OB Headrest (**STL-OTS250**) with Multi Axis Hardware (**HMO475**), Swing-Away ([**OTTOS9-L**] Left) with Fiber Optic Mount (**FO650**).

OB04

Sip & Puff Mount



Includes: Small OB Headrest (**STL-OTS250**) with Multi Axis Hardware (**HMO475**), Swing-Away ([**OTTOS9-L**] Left) with Sip and Puff Mount (**SP625**).