

Ultra Lightweight Pediatric Back

Standard Hardware Mount

The Ultra-Lightweight Pediatric True Forming Back (TFBN) is Stealth's Niño Seating Series latest design in comfort and custom fit options. This back system allows for optimal fitting through a flexible back structure. The lightweight, aluminum plate embraces the contour of the back while providing and maintaining the proper positioning of the trunk and pelvis.

Our hardware is durable and easily accessible. The growing kit allows for this new amazing back to have adjustability up to 2 inches in both height and width. The TFBN can be mounted using a one point center piece bracket, included with any package, or the traditional two point way with the use of a "Rigidizer Bar" when needed.

High quality stitch work guarantees long lasting resiliency. The TFBN comes standard with 1/2" medium SunMate® and 1/2" HR3 foam for ultimate comfort while maintaining the ultra-lightweight feature. Other combinations of foam are always available for optimal comfort, provided the maximum thickness does not exceed 1"



Single Slot Mount



Double Slot Mount

Additional Options



Rigidizer Bar

Ultra Lightweight Headrest And Hardware



Removable Cushion

Easy to install, washable and removable covers and foams. Standard 1/2" Medium SunMate and 1/2" HR3 Foam



We engineer custom mounting hardware

If it doesn't exist we customize it for you!



Ask about the ultimate true comfort the True Forming Back Series

Stealth Products
Beyond Expectations

FOLLOW US

www.stealthproducts.com
www.facebook.com/StealthProducts



1.800.965.9229

