

True Forming Back

break through
true forming back

The True Forming Back bestows absolute comfort, custom fit options, and trouble-free adjustability. Our hardware is durable, easily accessible, and is hidden with a Velcro attached back for a clean, sleek look.

The TFB allows for optimal fitting through a flexible back structure. The lightweight aluminum plate embraces the contour of the back providing and maintaining the proper positioning of the trunk and pelvis.

We offer high quality stitch work that guarantees long lasting resiliency. The standard model TFB comes with dual layer foam for ultimate comfort. Choosing from our assortment of types and range of densities in foam, each individual can be provided with the perfect comfort.



Single Wings



Double Wings



Lightweight Back
with integrated tracks for easy wings and lateral support installation



Available in a variety of sizes



J&L Mounting Hardware



ADI Mounting Hardware



Trouble Free Mounting Hardware



Foam In Place

The Foam-In-Place (FIP) kit offers the possibility to create a mold to increase the contact surface with the back of the client to correct posture.

*We engineer custom mounting hardware
If it doesn't exist we customize it for you!*



*Ask about the back that grows
Ultra Lightweight Pediatric Backs*

Stealth Products
Beyond Expectations

FOLLOW US

www.stealthproducts.com
www.facebook.com/StealthProducts



1.800.965.9229

True Forming Back (TFB)

PO#

Date: _____ Customer: _____ Location: _____
 Contact: _____ Email: _____
 Phone: _____ Fax: _____

Model Type: _____ Color: _____
 Frame Width: _____
 Cane Size: 3/4" 7/8" 1"

Contact Material: _____
 Non-Contact Material: _____
 Monogram: Block Script
 Name: _____
 Color: _____

Standard Sizes						
Frame	A	B	G	H	I	J
11" x 12"	7	13	8	2	3	7
13" x 14"	9	15	10.5	3	4	9.5
15" x 16"	11	17	11.5	3	4	10.5
17" x 18"	13	19	12.5	3	4	11
19" x 20"	15	21	13.5	3.5	4	12
21" x 22"	17	23	14.5	3.5	4	13
23" x 24"	19	25	15.5	3.5	4	14
25" x 26"	21	27	15.5	3.5	4	14
Custom						

Step 1

Foam In Place Kits

TFB100 4 Point Trouble Free Mounting Hardware with Drop Stops \$944.00

TFB101 4 Point Adjustable J Hook with Fixed/Rotating Hold Downs \$844.00

TFB102 4 Point Adjustable J Hook with Fixed Hold Downs \$844.00

TFB104 2 Point Quick Release Mounting Hardware \$1126.00
includes cover for Foam In Place and 1 Unit of FIPS Solution.

Standard Kits

TFB105 4 Point Adjustable J Hook with Fixed/Rotating Hold Downs \$759.00

TFB106 4 Point Adjustable J Hook with Fixed Hold Downs \$759.00

TFB108 4 Point Trouble Free Mounting Hardware with Drop Stops \$828.00

TFB109 2 Point Quick Release Mounting Hardware \$941.00
includes standard foam, 1/2" Medium SunMate over 1/2" HR3.

Step 2

Modifications

TFBM200 Extended Back: Width: Height: \$0.00

TFBM201 Split Wing And Cover \$325.00

BMD214 Aluminum Base For Back \$75.00

BMD227 Oversized, 22"-24" Either Width Or Height \$60.00

BMD228 Oversized, 25"-28" Either Width Or Height \$60.00

Step 3

Foam Options And Upgrades

[†]Standard foam is 1/2" Medium SunMate over 1/2".

1/2" SunMate M S XS Back: \$45.00 Lateral Pad: \$25.00

1" SunMate M S XS Back: \$75.00 Lateral Pad: \$35.00

1/2" T-Foam M S XS Back: \$55.00 Lateral Pad: \$25.00

1" T-Foam M S XS Back: \$85.00 Lateral Pad: \$35.00

BMD293 Spinal Relief \$75.00

Step 4

Additional Upgrades And Accessories

TWBL Lateral Hrdw. Replacement For Upper Wing On Split Wing \$225.00

TWBT Thin Lateral Hrdw. Replacement For Upper Wing On Split Wg \$225.00

LTP502 Lateral Thoracic Pad Flat Wood Width: Height: \$75.00



LTP506 Lateral Thc. Pad Curved Metal Width: Height: \$86.00

LTHW504 Lateral Thoracic All Season Mount \$95.00

LTHW511 Full Surface Contact Feature For TWBL Lateral Only \$75.00

LTHW511 Full Length Tracks \$33.00

Price shown for each unless notify



Order Reference Chart

Qty	Part #	Description	Price
Total			

