

True Forming Back (TFB)

PO#

Date: _____ Customer: _____ Location: _____
 Contact: _____ Email: _____
 Phone: _____ Fax: _____

Model Type: _____ Color: _____
 Frame Width: _____
 Cane Size: 3/4" 7/8" 1"

Contact Material: _____
 Non-Contact Material: _____
 Monogram: Block Script
 Name: _____
 Color: _____

Standard Sizes

| Frame | A | B | G | H | I | J |
|-----------|-----------|-----|-------|------|----|-------|
| 11" - 12" | 8" - 9" | 13" | 8" | 2" | 3" | 7" |
| 13" - 14" | 10" - 11" | 15" | 10.5" | 3" | 4" | 9.5" |
| 15" - 16" | 12" - 13" | 17" | 11.5" | 3" | 4" | 10.5" |
| 17" - 18" | 14" - 15" | 19" | 12.5" | 3" | 4" | 11" |
| 19" - 20" | 16" - 17" | 21" | 13.5" | 3.5" | 4" | 12" |
| 21" - 22" | 18" - 19" | 23" | 14.5" | 3.5" | 4" | 13" |
| 23" - 24" | 20" - 21" | 25" | 15.5" | 3.5" | 4" | 14" |
| 25" - 26" | 22" - 23" | 27" | 15.5" | 3.5" | 4" | 14" |

Step 1

Foam In Place Kits

- TFB100** 4 Point Trouble Free Mounting Hardware with Drop Stops \$944.00
- TFB101** 4 Point Adjustable J Hook with Fixed/Rotating Hold Downs \$844.00
- TFB102** 4 Point Adjustable J Hook with Fixed Hold Downs \$844.00
- TFB104** 2 Point Quick Release Mounting Hardware \$1126.00
includes cover for Foam In Place and 1 Unit of FIPS Solution.

Standard Kits

- TFB105** 4 Point Adjustable J Hook with Fixed/Rotating Hold Downs \$759.00
- TFB106** 4 Point Adjustable J Hook with Fixed Hold Downs \$759.00
- TFB108** 4 Point Trouble Free Mounting Hardware with Drop Stops \$828.00
- TFB109** 2 Point Quick Release Mounting Hardware \$941.00
includes standard foam, 1/2" Medium SunMate over 1/2" HR3.

Step 2

Modifications

- TFBM201** Split Wing And Cover \$325.00
- BMD301** Bent Aluminum Base w/ Built in Tracks For Back \$85.00
- BMD227** Oversized, 22"-24" Either Width Or Height \$60.00
- BMD228** Oversized, 25"-28" Either Width Or Height \$60.00

Step 3

Foam Options And Upgrades



[†]Standard foam is 1/2" Medium SunMate over 1/2."

- 1/2" SunMate** M S XS Back: \$45.00 Lateral Pad: \$25.00
- 1" SunMate** M S XS Back: \$75.00 Lateral Pad: \$35.00
- 1/2" T-Foam** M S XS Back: \$55.00 Lateral Pad: \$25.00
- 1" T-Foam** M S XS Back: \$85.00 Lateral Pad: \$35.00
- BMD293** Spinal Relief \$75.00

Step 4

Additional Upgrades And Accessories

- TWBL** Lateral Hrdw. Replacement For Upper Split Wing \$225.00
- TWBTl** Thin Lateral Hrdw. Replacement For Upper Split Wing \$225.00
- TWBSTL** Thin w/ Link Lateral Hrdw. Replacement For Upper Split Wing \$225.00
- LT520** S/A Modular Hrdw. Replacement For Upper Split Wing \$265.00
- LTP502** Lateral Thoracic Pad Flat Wood Width: Height: \$75.00
- LTP505** Lateral Thoracic Pad Flat Metal Width: Height: \$75.00
- LTP507** Lateral Thc. Pad Curved Wood Width: Height: \$86.00
- LTP506** Lateral Thc. Pad Curved Metal Width: Height: \$86.00
- LTHW509** Lateral Thoracic All Season Mount \$95.00
- LTHW511** Full Surface Contact Feature For TWBL Lateral Only \$75.00
- Tracks** Single Slot Adjustment Rail For Wood Base Only \$33.00
 - 2" Track 3" Track 4" Track 6" Track 8" Track
 - 10" Track 12" Track 14" Track 16" Track 18" Track



Aluminum Base



Aluminum Base with Split Wings

