


TECHNICAL SPECIFICATIONS



MTHW640

A hip abductor or pommel is used to maintain the legs in an abducted position. The MTHW640 fixed flip down hardware includes a fixed mounting bracket and swing

The MTHW640 includes:

1. Pad Mount
2. Height Adjustment
3. Flip Down Mechanism

DIMENSIONS



Maximum Hardware Height—4 3/4" Minimum Hardware Height—3 3/4"

Pad Mount Length—3" Hardware Length—4"

MEDIAL THIGH MODIFICATION



Adjustable Pad

Stealth's MTMD-42 allows adjustment in multiple angles/planes. Comes standard with dual 3x5" pads and slip cover.

TECHNICAL DATA

1. 5/32 tool for 1/4-20 button head screws (Mounting, height adjustment and pad mount screws).

SELLING POINTS

1. Flip down technology that makes user transfers simple and easy.
2. Height adjustable to better fit the seat and need of the user.

MEDIAL PADS



Wedge Style Pads

MTMD-12—Small Pad
MTMD-13—Medium Pad
MTMD-14—Large Pad
MTMD-19—Custom Pad



Oval Style Pads

MTMD-21—Extra Small Pad
MTMD-22—Small Pad
MTMD-23—Medium Pad
MTMD-24—Large Pad
MTMD-29—Custom Pad




Knee Cut Out Style Pads

MTMD-32—Small Pad
MTMD-33—Medium Pad
MTMD-34—Large Pad
MTMD-39—Custom Pad

Size	Length	Width	Height
Extra Small*	3 1/2"	2 1/2"	3"
Small	4"	3"	3"
Medium	5"	4 1/2"	3"
Large	5"	5"	4"


*Only oval style pads come in the extra small size.

WEDGE STYLE PAD DIMENSIONS




	Size	A	B	C	D
MTMD-12	Small	3"	4"	3"	2"
MTMD-13	Medium	4.5"	5"	3"	2"
MTMD-14	Large	5"	5"	4"	3"
MTMD-19	Custom	Custom Defined Dimensions			

OVAL STYLE PAD DIMENSIONS



	Size	A	B	C
MTMD-21	Extra Small	2.5"	3.5"	3"
MTMD-22	Small	3"	4"	3"
MTMD-23	Medium	4.5"	5"	3"
MTMD-24	Large	5"	5"	4"
MTMD-29	Custom	Custom Defined Dimensions		

KNEE CUT-OUT PAD DIMENSIONS



	Size	A	B	C	D
MTMD-32	Small	3"	4"	3"	2"
MTMD-33	Medium	4.5"	5"	3"	2"
MTMD-34	Large	5"	5"	4"	3"
MTMD-39	Custom	Custom Defined Dimensions			

ADJUSTABLE POMMEL MODIFICATION SPECIFICATIONS

This adjustment screw allows you to change the angle of the pads.

The pad connects to the planar mount with a track and track nut. This allows the pad to be positioned fore and back on the planar mount.

These screws allow the hardware's depth position to be adjusted.

The MTMD-42 mounts onto a medial thigh support hardware with this bracket.



MTMD-42

Stealth's adjustable medial pad allows adjustment in multiple angles/planes. The Purpose of this adjustment is to allow for a more natural fit and comfort for the client. This modification is compatible with the MTHW-600, MTHW-610, MTHW-620, MTHW-630, and MTHW-640; and comes standard with dual 3x5" pads and slip cover.

ADJUSTABILITY



SLIP COVER

