



Where Style Meets Tech

i-Drive is the most advanced head array ever built, not only because of the wide positioning options offered by Stealth Products, but by the technology engineered behind it!



More Information
1 (800) 965-9229

New Technology for a New Mobility Era

Integrated circuit and digital signal processing has advanced significantly since the invention of the modern microprocessor. Today's smart phones and tablet devices are driving the technology in this area, making it affordable to apply sophisticated technologies to the mobility industry. The technology in the i-Drive product line makes use of these technologies by combining reliable integrated circuits with modern printed circuits board design and signal processing.



Easy Set Up App

The i-Drive's advanced programming software is not required to operate the drive control, but is available to ensure a customized fit and driving experience. Programming the i-Drive is simple. Simply connect a Windows PC, Laptop or Tablet and use the application for customization.



Like Us! :)



Custom Settings

With six available input ports, the i-Drive interface gives access to four proximity sensors and four mechanical switches allowing for eight independent function switches.

i Position

Hundreds of positioning options from our vast catalog...

i Set Up


The i-Drive App allows you to set your experience...

i Drive

Empowering you to accomplish your dreams...

Unique Features:

- **Dynamic Channel Assignment** - The "i-Drive" not only has the ability to manually plug and unplug from it's interface, it also have the ability to assign the port you plug into, any port can control any of the 5 outputs.
- **Double Tap Timing** - This is a unique feature that allows you to choose the speed of a "mode select command". Extends chair double tap limits.
- **Mode Selection** - Can be configured on any port and allows for unique placement of reset/mode changes.
- **Sensor Engagement Setting** - Allows for the highest level of driving experience for the client. Each sensor can be independently adjusted for sensitivity or activation control.
- **Bluetooth Technology** - (Soon) You will be able to program and control the mobility base via Bluetooth.



www.stealthproducts.com/idrive

Stealth Products

Powered by
TRIDENT
RESEARCH
Our Nation trusts them, So does Stealth!